

US (ultrasilicon)

NEW

MASSAGINGFIBRE

Every movement triggers a micromassage providing well being and lightness. Patented fibre.



CHOCABSORBER

The 360° compression around the whole calf absorbs the "shock wave" effect on muscles, veins and joints.



ULTRALIGHT

Wet or dry, COMPRESSPORT® probably remains the most lightweight (does not retain water).



ULTRARESIST

In normal use COMPRESSPORT® is tear-proof.



QUALITY GUARANTEE

COMPRESSPORT® guarantees the exchange of your product. (see conditions on our website <http://www.compressport.com>)



DTOX



Studies show a significant reduction in lactic acid during effort.

ANTI-FATIGUE



The shock and vibration absorption on muscles, veins and joints increases resistance to effort.

RECOVERY



No pain after strain! COMPRESSPORT® is used for recovery after showering (in a plane, at the office, at home)!

NEW

SILICONABSORBER

US® is the latest in the COMPRESSPORT® range. With its 220 soft silicon spikes which fully surround the ankle over a height of 7cm, US® optimises shock absorption and guarantees ideal support to the Achilles tendon throughout strenuous effort (helping absorb pain and tendon problems). The COMPRESSPORT® US® is a major leap forward in the compression world.


DESIGN & QUALITY SWISS COMPANY


2 SWISS COMPANY GARANTEE

By purchasing your COMPRESSPORT® you receive a warranty. You are required to validate your warranty on our website www.compressport.com category security. Registration is required for validation of your warranty. In return you receive a numbered certificate of guarantee that you can easily maintain and use if needed. Sportingly Team Compressport.

30 MONEY-BACK GUARANTEE

Test your COMPRESSPORT® for 30 days! Discover the performances of your COMPRESSPORT® at the track, on road, on trail or bike. Make your own opinion during 30 days. Also discover RECOVERY after effort with COMPRESSPORT®. Our guarantee, we certify to reimburse you, if you are not satisfied.

ADVANTAGES +

The Compressport US® is your "Recovery" partner after intense effort, allowing you to be fresh for your next training session. Virtually indestructible, your Compressport US® can be stretched several times its size and recover its original shape. Compressport US® resists both chlorine, or acidity due to perspiration and salt. Use your Compressport US® "incognito" with long pants or "fashion" as the "pros" in shorts and thongs. Being footless your Compressport US® is more hygienic (no need to wash each session), which prolongs its life.